

Sodium Acrylate 38% (w/w) – Preparation Protocol

Purpose

To prepare ~180 mL of 38% (w/w) sodium acrylate solution by neutralizing acrylic acid with sodium hydroxide under controlled temperature.

Safety

WARNING: This procedure uses **concentrated acrylic acid** (corrosive, lachrymator) and **10 N NaOH** (highly caustic). Both can cause severe burns and permanent eye damage within seconds.

- Perform **all steps in a certified fume hood**.
 - Wear **lab coat, safety goggles/face shield, and appropriate chemical-resistant gloves**.
 - Use a **spill kit** and have **eyewash and safety shower** accessible.
 - Neutralize and dispose of waste according to institutional regulations.
-

Reagents

- Acrylic Acid (Sigma-Aldrich, Cat. No. 8.00181)
- 10 N NaOH
- Deionized or double-distilled water (ddH₂O)

Equipment & Materials

- Glass beaker (250 mL or larger)
 - Glass bottle (≥100 mL) for NaOH
 - Magnetic stir plate and stir bar
 - Peristaltic pump with compatible tubing
 - pH meter with temperature probe
 - Ice bucket / thin plastic container for ice bath
 - Centrifuge (for later clarification), capable of 1000 × g
 - Storage bottles (screw-cap, chemical-resistant)
-

Target Preparation

- Final volume: **≈180 mL**
 - Composition (approximate):
 - 50 mL acrylic acid
 - 75 mL 10 N NaOH
 - 58 mL ddH₂O
-

Procedure

1. Initial Mixing

1. In a clean glass beaker, add **58 mL ddH₂O**.
 2. **Slowly** add **50 mL acrylic acid** to the water **while stirring**.
 - **Important:** Always add *acid to water*, not water to acid.
 3. Add a magnetic stir bar and stir for ~1 min to obtain a homogeneous solution.
-

2. Ice Bath & pH Probe Setup (*see setup picture*)

1. Place ice into a thin plastic container.
 2. Insert the beaker (water + acrylic acid) into the ice, then add water around it to form an ice–water bath.
 - The beaker should be **approximately half immersed** in the ice–water mixture.
 - Add more ice as needed to keep the bath cold and stable.
 3. Place the beaker on a magnetic stirrer and **start stirring**.
 4. Install the pH/temperature probe in the solution, ensuring:
 - The probe is **fully submerged** in the liquid.
 - It does **not interfere** with the stir bar.
-

3. NaOH & Peristaltic Pump Setup (*see setup picture*)

1. In a separate glass bottle, **pour 75 mL of 10 N NaOH**.
 2. Connect the peristaltic pump tubing:
 - Place the **inlet/collector end** at the **bottom** of the NaOH bottle.
 - Secure the **outlet end** above the beaker in the ice bath.
 3. Position the outlet so that it drips **into the beaker** but:
 - The tube **does not touch** the liquid.
 - There is enough clearance for the volume to increase during addition.
-

4. Controlled Neutralization

1. Start the peristaltic pump at **~50% of its maximum rate** to begin NaOH addition.
 - The solution temperature will start to rise.
2. When the solution temperature reaches **25 °C**, **stop the pump**.
3. Restart the pump at **~10%** and closely monitor **both temperature and pH**.
4. Adjust the pump speed (typically **around 14%**) so that:
 - The temperature stays **near 25 °C** (avoid >30 °C).
 - The NaOH is added **gradually**, preventing overheating and overshooting the pH.
5. As the pH approaches **7.0**, **slow the pump further** and then:
 - **Stop the pump immediately** once the pH is slightly above neutral (pH 7.0–7.5).

If pH Overshoots (e.g., pH \approx 10.5)

1. Remove the beaker from the ice bath and place it on the stir plate at **room temperature**.
2. Allow the temperature to stabilize while stirring.
3. Carefully add **small aliquots of acrylic acid** (e.g., 0.5–1 mL at a time) while monitoring pH, until it drops to **pH 7–8**.
4. Allow the solution to mix thoroughly after each addition before re-measuring pH.

5. Equilibration, Clarification & Storage

1. Transfer the neutralized solution to a **closed container** (loosely capped if still slightly warm) and allow it to stand **overnight at room temperature** or 4 °C to equilibrate.
2. The next day, inspect the solution:
 - If **no precipitate** is visible, proceed to storage.
 - If **precipitate is present**, centrifuge at **1000 × g for 10 min** and collect the **clear supernatant**.
3. Aliquot the clear sodium acrylate solution into storage bottles.
4. **Storage recommendations:**
 - **Preferred:** –20 °C for long-term storage.
 - **Alternative:** 4 °C for shorter term use.
5. Label each container clearly with:
 - “Sodium acrylate 38%”
 - Date of preparation
 - Preparer’s initials
 - Storage conditions and any relevant hazard warnings.

Notes & Troubleshooting

- **pH target:** A final pH of **7.0–8.0** is typically suitable for downstream applications; adjust if your protocol requires a specific value.
- **Foaming:** If foaming occurs, reduce stirring speed while maintaining good mixing.
- **Ice bath:** If temperature creeps up, add more ice and/or reduce the pump speed.

Picture of the set up:

